

Cytiva Xcellerex™ XDR-2000 - approximately 400-2000 L working volume

Local source PDF | preset cytiva-xdr-2000 | compiled 10 September 2026

This PDF is stored inside the GitHub project and records the evidence used by the simulator. It is not a copied vendor brochure or manual. The official manufacturer references below remain the authoritative originals.

Manufacturer: Cytiva

Product-contact material: Plastic / single-use

Working volume: 400 L to 2000 L

Nominal vessel or bag size: 2000 L

Simulation profile: mammalian

Evidence status: Official family evidence with derived minimum working volume

Original manufacturer references

1. Cytiva - Xcellerex XDR bioreactor system catalogue

[Open official manufacturer source online](#)

<https://www.cytivalifesciences.co.jp/catalog/1754.html>

Supports: XDR family, XDR-10 vessel parameters and larger-scale model range.

2. Cytiva - XDR family scale information

[Open official manufacturer source online](#)

<https://info.cytivalifesciences.com/bpinsights2019.html>

Supports: Supporting family statement covering the XDR scale range.

Evidence basis

Volume: Nominal XDR model sizes are source-supported. Minimum working volumes for XDR-50 and larger are calculated at approximately 20% of maximum scale.

Material: Disposable single-use bag is the product-contact vessel.

- Large XDR presets use a generic mammalian stirred-tank profile unless a model-specific field is stored.

Record note

This source PDF is paired with the detailed HTML provenance record in the same repository. Simulator-specific defaults and derived assumptions are listed there.

Copyright and trademarks in the linked manufacturer documentation remain with the respective vendor.