

Mobius(R) iFlex 1000 L Bioreactor - approximately 200-1000 L working volume

Local source PDF | preset merck-iflex-1000 | compiled 10 September 2026

This PDF is stored inside the GitHub project and records the evidence used by the simulator. It is not a copied vendor brochure or manual. The official manufacturer references below remain the authoritative originals.

Manufacturer: Merck

Product-contact material: Plastic / single-use

Working volume: 200 L to 1000 L

Nominal vessel or bag size: 1000 L

Simulation profile: mammalian

Evidence status: Family evidence with inherited/derived working range

Original manufacturer references

1. Merck - Mobius single-use bioreactors / iFlex family

[Open official manufacturer source online](#)

<https://www.merckmillipore.com/LT/en/products/pharma-and-biopharma-manufacturing/bioprocessing-cell-culture/single-use-bioreactors>

Supports: Current Mobius portfolio, iFlex family scale and single-use architecture.

2. Merck - Mobius single-use bioreactor scale announcement

[Open official manufacturer source online](#)

https://www.merckmillipore.com/AR/en/20131223_191341

Supports: Supporting historical working-volume ranges for 50, 200, 1000 and 2000 L Mobius scales.

Evidence basis

Volume: The iFlex family and nominal scale are source-supported. The working-volume range is based on official Mobius scale information and the expected family turndown; confirm against the current configuration data sheet.

Material: Disposable polymer bag is the product-contact vessel; reusable support hardware is not classified as the vessel material.

- These presets are primarily mammalian-cell-culture archetypes.

Record note

This source PDF is paired with the detailed HTML provenance record in the same repository. Simulator-specific defaults and derived assumptions are listed there.

Copyright and trademarks in the linked manufacturer documentation remain with the respective vendor.