

Biostat(R) RM with Flexsafe(R) RM 1 L bag - 0.1-0.5 L working volume

Local source PDF | preset sartorius-rm-1 | compiled 10 September 2026

This PDF is stored inside the GitHub project and records the evidence used by the simulator. It is not a copied vendor brochure or manual. The official manufacturer references below remain the authoritative originals.

Manufacturer: Sartorius

Product-contact material: Plastic / single-use

Working volume: 100 mL to 500 mL

Nominal vessel or bag size: 1 L

Simulation profile: rocking

Evidence status: Direct official product-family evidence

Original manufacturer references

1. Sartorius - Biostat RM and Flexsafe RM Bags

[Open official manufacturer source online](#)

<https://www.sartorius.com/en/products/fermentation-bioreactors/single-use-bioreactors/biostat-rm-flexsafe-rm>

Supports: Bag sizes, working-volume ranges and rocking-bag architecture.

Evidence basis

Volume: Bag size and working-volume range follow the official Biostat RM / Flexsafe RM matrix.

Material: The product-contact vessel is a disposable rocking bag; the hardware support is not the wetted vessel.

- The simulator uses a rocking profile. Impeller-related fields are computational placeholders and not a claim that the bag contains an impeller.

Record note

This source PDF is paired with the detailed HTML provenance record in the same repository. Simulator-specific defaults and derived assumptions are listed there.

Copyright and trademarks in the linked manufacturer documentation remain with the respective vendor.