

Thermo Scientific™ HyPerforma™ 5:1 S.U.B. 1000 L - 200-1000 L working volume

Local source PDF | preset thermo-sub-1000 | compiled 10 September 2026

This PDF is stored inside the GitHub project and records the evidence used by the simulator. It is not a copied vendor brochure or manual. The official manufacturer references below remain the authoritative originals.

Manufacturer: Thermo Fisher Scientific

Product-contact material: Plastic / single-use

Working volume: 200 L to 1000 L

Nominal vessel or bag size: 1000 L

Simulation profile: mammalian

Evidence status: Official family evidence with derived 5:1 minimum

Original manufacturer references

1. Thermo Fisher Scientific - Single-use bioreactor portfolio

[Open official manufacturer source online](#)

<https://www.thermofisher.com/es/es/home/life-science/bioproduction/single-use-bioprocessing/single-use-equipment/single-use-bioreactors.html>

Supports: Supporting product-family overview and available HyPerforma single-use scales.

2. Thermo Fisher Scientific - HyPerforma 5:1 S.U.B. 50 L

[Open official manufacturer source online](#)

<https://www.thermofisher.com/order/catalog/product/es/es/SUB05008400>

Supports: 50 L model and 5:1 single-use bioreactor platform.

Evidence basis

Volume: Nominal model sizes and the 5:1 platform are source-supported. The minimum working volume is calculated as nominal volume divided by five.

Material: Disposable BPC is the product-contact vessel; stainless-steel support hardware is not the wetted vessel.

- These presets use a generic mammalian-cell-culture profile.

Record note

This source PDF is paired with the detailed HTML provenance record in the same repository. Simulator-specific defaults and derived assumptions are listed there.

Copyright and trademarks in the linked manufacturer documentation remain with the respective vendor.