

Thermo Scientific™ HyPerforma™ S.U.F. 30 L - 6-30 L working volume

Local source PDF | preset thermo-suf-30 | compiled 10 September 2026

This PDF is stored inside the GitHub project and records the evidence used by the simulator. It is not a copied vendor brochure or manual. The official manufacturer references below remain the authoritative originals.

Manufacturer: Thermo Fisher Scientific

Product-contact material: Plastic / single-use

Working volume: 6 L to 30 L

Nominal vessel or bag size: 30 L

Simulation profile: microbial

Evidence status: Direct official model evidence

Original manufacturer references

1. Thermo Fisher Scientific - HyPerforma S.U.F. 30 L

[Open official manufacturer source online](#)

<https://www.thermofisher.com/order/catalog/product/SUF0030.9003>

Supports: 6-30 L working range, vessel geometry, impeller configuration and agitation range.

Evidence basis

Volume: The 6-30 L working-volume range is directly supported.

Material: Disposable BPC is the product-contact vessel inside reusable support hardware.

- The 3:1 reactor H/D, three Rushton-style impellers and 600 rpm maximum are model-specific source-supported fields.

Record note

This source PDF is paired with the detailed HTML provenance record in the same repository. Simulator-specific defaults and derived assumptions are listed there.

Copyright and trademarks in the linked manufacturer documentation remain with the respective vendor.